




CANDELA®



DERMA LASER
Your Laser Service Partner

Classic V-BEAM

VBEAM Specifications



Laser Type	Pulsed Dye
Wavelength	595 nm
Repetition Rate	1.5 Hx maximum
Pulse Duration	0.45-40 ms (user-selectable)
Spot Sizes	5 mm—up to 15 J/cm
(with	7 mm—up to 15 J/cm
maximum	10 mm—up to 7.5 J/cm
energy)	3 x 10 mm—up to 25 J/cm
Dynamic Cooling	Integrated
Device	Spray: 20-100 ms
	Delay: 10-100 ms
	High Volume Canister (1,000 g)
Unit Colors	Available in 4 choices:
	Soft Rose 
	Ice Blue 
	Vintage Violet 
	Cool Gray 
Dimensions	43.5"H x 18"W x 32"D
	(111 cm x 46 cm x 81 cm)
Weight	300 lbs (136 kg)
Electrical	230 VAC 50/60 Hz, 20 A


TREAT VIRTUALLY ALL VASCULAR LESIONS

Vbeam's user selectable pulse durations and spot sizes are the key to its versatility. With the addition of new 5 mm spot size and a short 0.45 ms pulse duration, you can treat a greater range of vascular conditions than ever before, without patient downtime caused by purpura.

- Rosacea
- Facial telangiectasia
- Leg telangiectasia
- Periorbital wrinkles
- Angiomas
- Hemangiomas
- Port wine stains
- Poikiloderma of Civatte
- Scars
- Stretch marks
- Warts
- Venous lakes
- Psoriasis
- Benign gynecological vascular lesions
- Other benign cutaneous vascular lesions

For Sale & Service Contact 0321-9455593

PRE AND POST TREATMENT RESULTS



Brown Pigmentation



Brown Pigmentation




Periorbital Wrinkles



Leg Veins




Leg Veins



Port Wine Stain




Hemangioma



Nose Veins




Poikiloderma of Civatte



Acne Redness



Rosacea



Rosacea



Port Wine Stain



Venous Lake



Scar



Warts