



PERFECT FOR EVERYTHING PERFECTA

Perfecta can do it all: from rejuvenation by correcting red and brown skin discoloration, like Aesthetica, to eliminating virtually all vascular and pigmented lesions, like the Platinum. The Perfecta features advanced micro-pulse technology, multiple handpieces for treating various spot sizes, and our patented Dynamic Cooling Device™ (DCD™) cooling system, which all but eliminates patient downtime by preventing purpura, or bruising. And like all Candela products, Perfecta is backed by Candela's unmatched training, marketing and service—including a financing program that's right for your practice.

One laser for all these treatment conditions:

- Skin Rejuvenation
- Diffuse Redness
- Pigmentation
- Freckles
- Sun Spots
- Rosacea
- Acne
- Facial Veins
- Leg Veins
- Pigmented Lesions
- Wrinkles
- Poikiloderma of Civatte
- Angiomas
- Hemangiomas
- Port Wine Stains
- Scars and Stretch Marks
- Warts
- Venous Lakes
- Psoriasis
- Benign Gynecological Vascular Lesions
- Other Benign Cutaneous Vascular Lesions



Vbeam PERFECTA SPECIFICATION

Laser Type	Pulsed Dye
Wavelength	595 nm
Pulse	0.45-40 ms
Spot Sizes (diameter) with maximum energy	3 mm - 40 J/cm ² 5 mm - 30 J/cm ² 7 mm - 20 J/cm ² 10 mm - 10 J/cm ² 12 mm - 7 J/cm ² 3x10 mm - 25 J/cm ² 7 mm PL - 15 J/cm ² 10 mm PL - 10 J/cm ²

Beam Delivery	Lens-coupled 3 m optical fiber with handpiece
Pulse Control	Fingerswitch, Footswitch
Dimensions	43" H x 16" W x 30" D (109 cm x 40 cm x 76 cm)
Weight	290 lbs (132 kg)
Electrical	230 V 50/60 Hz, single phase, 18A

For Sale & Service Contact 0321-9455593

PERFECT FOR EVERYTHING



PRE AND POST TREATMENT RESULTS



Brown Pigmentation



Brown Pigmentation




Periorbital Wrinkles



Leg Veins



Leg Veins



Port Wine Stain




Hemangioma



Nose Veins




Poikiloderma of Civatte



Acne Redness




Rosacea



Rosacea



Port Wine Stain



Venous Lake



Scar



Warts

