

# Only Smartlipo TriPlex offers 3 wavelengths

Including 1440 nm for faster disruption of adipose tissue – delivered by multiplex sequencing and intelligent delivery systems. It's laser lipolysis redefined, only from Cynosure. Why settle for anything less?

**Greater absorption.**

1440 nm has 40 times greater absorption in fatty tissue than diode lasers.

**More choice.**

Multiplexed combination wavelength therapy for optimal fat disruption, high-definition sculpting capability and tissue tightening through coagulation.

**Better Results.**






SmartSense, ThermoGlide and ThermoView intelligent delivery systems help you treat with the maximum dose, safely.



Smartlipo TriPlex™  
Aesthetic Workstation

System Type	Nd:YAG		
Wavelength	1064nm	1320nm	1440nm
Pulse Width (Max)	500µs	1000µs	1000µs
Power Output Options	40W	24W	15W
Max Rep Rate	40Hz	40Hz	25Hz
	MultiPlexing Options		
	1064nm and 1320nm	1064nm and 1440nm	
	500µs and 1000µs	500µs and 1000µs	
	40Hz	40 Hz and 25Hz	
Aming Beam	Red		
Electrical Requirements	230 VAC/16A		
Size (HxWxD)	37.4" x 19.7" x 25.6" (95cm x 50cm x 65cm)		
Weight	220 lbs (100kg)		

For Sale & Services Contact: 0321-9455593



Dramatic  
Results in  
Less Time.

No wonder  
Cynosure's Smartlipo  
multiplexing  
technology is the  
leading choice for  
laser lipolysis around  
the world.